



Title: NORC Portal User Guide

MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	1 of 50

NORC Portal User Guide

Clinical Site Guide to the NORC Site Portal

How to Use This Guide

This guide is written for clinical site staff who use the Netramind ophthalmic reading center's (NORC) Site Portal every day. You can read the chapters in order when you are first getting started, or use this document as a reference to assist with your workflow

Portal features and requirements may vary by study. Refer to the study protocol and applicable study materials to confirm which features are enabled.



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	2 of 50

Title: NORC Portal User Guide

Table of Contents

How to Use This Guide 1

Table of Contents..... 2

1. Introduction 5

1.1 Purpose..... 5

1.2 Intended Audience 5

1.3 Site and Reading Center Responsibilities 5

1.4 Prerequisites..... 5

2. Getting Started 6

2.1 Create a Site User Account (Sign Up) 6

2.2 Sign In (Login + MFA) 7

2.3 Reset Password 8

2.4 Sign out..... 9

2.5 Session timeout (Idle Warning and Auto Logout) 10

3. Portal Navigation 10

3.1 Assigned Studies..... 10

3.2 Header and Sidebar 11

3.3. Dashboard 11

3.3.1 What Needs Attention Panel 12

4. Site Readiness 14

4.1. Study documents..... 14

4.2. Site certification 15

4.2.1 Certification Statuses 15

4.2.2 Recommended Certification Pathway 17

4.2.3 Site and Device Information..... 17

4.2.4 Test Upload 18

4.2.5 Attestation and Submission of Images..... 19

5. Subjects..... 20

5.1 Subjects List 20

5.2 Register a Subject..... 21



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	3 of 50

Title: NORC Portal User Guide

5.3 Subject Details and Status Updates	22
5.4 Subject Visits	23
6. Imaging submission	24
6.1 Prerequisites checklist.....	24
6.2 New Submission Wizard.....	25
6.3 Upload Imaging Files	26
6.4 Partial Submissions**.....	27
6.5 Review, Attest, and Submit Imaging	27
6.6 Replacement Images	29
6.7. Upload History and Submission Details	29
7. Automatic Slot Check.....	30
7.1 What you will see on imaging upload	31
7.2 Failure reasons you may see	32
7.3 Submitting while Slot Check is still running or has failed	32
8. Automatic Grade Check.....	33
8.1 What you will see on imaging upload	34
8.2 Submitting while Grade Check is still running or has failed.....	35
8.3 Grade Check tab (left menu)	35
9. Queries.....	36
9.1 Queries List.....	36
9.2 Respond to Queries.....	37
10. Data Dumps.....	38
10.1 Create and Upload a Data Dump	38
11. Support and Communications.....	40
11.1 Study Contacts.....	40
11.2 Help (FAQs and Contact support)	40
11.3. Notifications	41
11.4 Email notifications.....	42
12. Revision History.....	44
13. Signature Page.....	45
Appendix A — Key Status reference	46



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	4 of 50

Title: NORC Portal User Guide

A.1 Site Certification	46
A.2 Imaging Submissions	46
A.3 Queries	47
Appendix B — Glossary	47
Appendix C — Support and troubleshooting	48
C.1 Support contacts	48
C.2 Common issues	49
C.3 Related study materials to reference.....	50



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	5 of 50

Title: NORC Portal User Guide

1. Introduction

1.1 Purpose

This guide explains how clinical sites use the **NORC Site Portal** to complete site certification, register subjects, submit ophthalmic imaging data, respond to Reading Center queries, and—when the study allows it—transfer bulk data packages.

1.2 Intended Audience

Site investigators, coordinators, photographers, and other clinical site staff who have been assigned to a study or assist with imaging upload/query resolution in the NORC Site Portal.

1.3 Site and Reading Center Responsibilities

Clinical Site Users	Reading Center (NORC)
Create an account and sign in	Approves your account
Complete site certification (when required)	Reviews and certifies the site
Acknowledge study documents (when required)	Publishes documents and contacts
Register subjects and submit imaging	Performs technical QC on submissions
Respond to queries	Opens and closes queries
Submit data dumps (when enabled)	Receives and processes transfers

1.4 Prerequisites

Before you begin, make sure you have:

- ☐ A workstation with a supported modern web browser and a reliable internet connection
- ☐ Access to the email inbox you will use for portal verification codes and study notifications
- ☐ The study imaging requirements (Imaging Acquisition Manual and related documents), as provided for each study



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	6 of 50

Title: NORC Portal User Guide

2. Getting Started

2.1 Create a Site User Account (Sign Up)

Purpose: Request access to the NORC Site Portal. New accounts will remain pending until an administrator from the Reading Center approves them.

Steps:

1. Open the portal norc.netramind.ai. (Figure 1)
2. Select **Create account**. (Figure 2)
3. Enter your email address and request a verification code.
4. Enter the one-time code from your email.
5. Complete your profile: name, password, security questions, and role.
6. For a clinical site account, enter your **site name** and **location** when prompted.
7. Submit the registration. You will see a message confirming that your account is awaiting approval.




Figure 1. NORC portal landing page



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	7 of 50

Title: NORC Portal User Guide




Figure 2. Create account page

Notes:

- ☐ Choose the clinical site option (not sponsor or grader) when signing up as site staff. This will prevent any delays in reviewing your access request.
- ☐ You cannot sign in until the Reading Center sets your account to **active**. You should receive an email when your access is reviewed by an administrator.
- ☐ If approval is delayed, contact your Reading Center project contact (see [Appendix C](#)).

2.2 Sign In (Login + MFA)

Purpose: Securely access your account using your email and password, followed by a multi-factor authentication process involving a one-time verification code sent to your registered email.

Steps:

1. Open the portal login page. (Figure 3)



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	8 of 50

Title: NORC Portal User Guide

2. Enter your registered email and password.
3. Enter the verification code from your email.
4. After successful multi-factor authentication (MFA), you land on **Assigned Studies**.




Figure 3. Login/sign-up MFA page

Notes:

- ☐ Ensure you have access to the registered email address linked to your account, as verification codes are sent to this address and remain valid for 10 minutes.
- ☐ If login fails after your account has been approved, confirm that you are using the approved email address and contact support if the problem persists.

2.3 Reset Password

Purpose: Regain access to your account if you forget your password. Password reset is available for **active** accounts only.

Steps:



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	9 of 50

Title: NORC Portal User Guide

1. On the login page, select **Forgot password**.
2. Enter your email and request a verification code.(Figure 5)
3. Enter the code from your email. (Figure 4)
4. Answer your security questions.
5. Set a new password and confirm it.
6. Sign in with the new password (including MFA).




Figure 4. One-Time code email




Figure 5. Reset password

2.4 Sign out

Open the user menu in the header (your name or avatar area) and select **Sign Out**.



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	10 of 50

Title: NORC Portal User Guide

2.5 Session timeout (Idle Warning and Auto Logout)

After a period of inactivity, the portal displays a **Still there?** warning with a countdown. Select **Stay signed in** to continue your session. If no action is taken before the countdown expires, you will be signed out automatically and returned to the login page.

Notes:

- It is generally recommended to finish uploads or other multi-step work before stepping away for an extended period.

3. Portal Navigation

3.1 Assigned Studies

After signing in, **Assigned Studies** displays the studies available to your site. Select **Open Study** to access the applicable study. (Figure 6)




Figure 6. Assigned studies page



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	11 of 50

Title: NORC Portal User Guide

3.2 Header and Sidebar

Inside the main page of a study (Figure 7), the **NORC Site Portal** header displays:

- ☐ Brand / home link back toward studies
- ☐ **Study switcher** (when you have more than one study)
- ☐ Your **site name** and **site number**
- ☐ **Notifications** (bell icon)
- ☐ **User menu** (including Sign Out)

The left sidebar includes:

Menu item	Use
Dashboard	Study overview, summary counts, and actionable attention items
Certification	Site certification package
Subjects	Subject list and registration
New Submission	Start an imaging submission
Upload History	Past and in-progress submissions
Data Dumps	Bulk zip transfers (only if enabled for the study)
Queries	Reading Center queries
Grade Check	Track grade check status for imaging submissions
Documents	Study documents and acknowledgment
Study Contacts	Reading Center and support contacts directory
Help	In-app FAQs and Contact support message form

3.3. Dashboard

Purpose: The **Dashboard** provides an overview of your site's activity and study status, including active subjects, upcoming or overdue uploads, submissions pending QC, and open queries. Select a summary card to view the related items. (Figure 7)

Depending on your site and study status, the Dashboard may also provide quick access to **Site Certification**, **Study Documents**, **Subjects**, **New Submissions**, and **Data Dumps** (when enabled). (Figure 9)



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	12 of 50

Title: NORC Portal User Guide




Figure 7. Study dashboard

3.3.1 What Needs Attention Panel

Purpose: The **What needs attention** section is an actionable work queue. Each row is a link that takes you directly to the task—unlike the summary cards, which display item counts.

Typical attention items include:

Attention item (examples)	What to do
Respond to query ...	Open the query and submit your response
... upload overdue — ...	Start or complete the overdue visit upload
... upload due soon — ...	Start the visit upload that is approaching its deadline
Finish ... submission — ...	Resume a Draft submission or continue a Partially Submitted visit

The list shows a list of action items ordered by highest-priority items that need attention. (Figure 8)



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	13 of 50

Title: NORC Portal User Guide

 **What needs attention** · 14

Respond to query Q-003 — Month 12 (Test301-101-0001)

Month 8 upload overdue — Test301-101-0004

Screening upload due soon — Test301-101-0007

Screening upload due soon — Test301-101-0008

Screening upload due soon — Test301-101-0010

Grade check failed — Optional Imaging (OD) · Month 4 (Test301-101-0001)

Grade check failed — Optical Coherence Tomography (OD) (6x6 mm vol scan, 97 B-Scans) · Screening (Test301-101-0001)

[View all](#)

Figure 8. What needs attention card on study dashboard

Studies > Test-301 > Dashboard

Test-301 Dashboard

Test - GA Study Phase 3

Register Subject

Sponsor
Test Pharma

Site
Mercy (101)

Activation / Certification
Active Not Applicable

Imaging manual / Contact
v1.0
Julie Tran (NORC - Co-ordinator)

ACTIVE SUBJECTS
10

VISITS AWAITING UPLOAD
4

UPLOADS PENDING QC
7

OPEN QUERIES
1

LATE SUBMISSIONS
1

Figure 9. Data cards and site certification status on study dashboard

Notes:

THIS DOCUMENT CONTAINS CONFIDENTIAL AND PROPRIETARY INFORMATION – UNCONTROLLED COPY IF PRINTED



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	14 of 50

Title: NORC Portal User Guide

- ☐ Open queries and overdue or due-soon visits are good daily checks during active enrollment and imaging windows to ensure resolution in a timely manner.
- ☐ Certification, document acknowledgment, and site activation are handled through the **header buttons** and **amber banners**, not through the “What Needs Attention” list.
- ☐ The companion **Recent activity** panel summarizes recent events at your site. When certified, you can also open **View Upload History** from that area.

4. Site Readiness

Purpose: Complete readiness work before any study subject registration and production imaging uploads. Exact requirements are **per study design**.

Sign in

- Open study
- Complete Certification (if required)
- Acknowledge Documents (if required)
- Wait until site status is Active (Reading Center)
- Register subjects
- Submit imaging / respond to queries

Sites and associated users cannot register subjects or upload imaging until the site is marked as **Active** in the portal, which may require additional certifications if required per protocol. Readiness banners on the Dashboard also provide reminders when activation, certification, or document acknowledgment is still outstanding for your site.

4.1. Study documents

Purpose: Access study-related documents, including imaging manuals and reference PDFs. When acknowledgment is **Required** for the study, you must acknowledge documents before production uploads. When policy is **Not Applicable**, acknowledgment does not gate uploads (**as per study design**).

Steps:

1. Open **Documents**. (Figure 10)
2. Open each listed document and review it.
3. If acknowledgment is required, confirm acknowledgment (including name and password when prompted).
4. Confirm that documents show as **Acknowledged**.



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	15 of 50

Title: NORC Portal User Guide



Figure 10. Documents tab showing available documents for selected study

Notes:

- ☐ The latest Imaging Acquisition Manual version aligned with a specific study will be listed under Documents.
- ☐ If uploads are blocked, check that required documents are acknowledged and that the site is **Active**.
- ☐ When a new study document is published that requires acknowledgment, your site may receive an in-app notification.

4.2. Site certification

Purpose: Confirm site readiness for production imaging by completing the required personnel, devices, test images, and attestation requirements. Certification **requirements are study-specific and may be marked Not Applicable** when not required. (Figure 11, Figure 12)

4.2.1 Certification Statuses

Certification status	Meaning
Not Started / In Progress	Complete the steps on screen
Pending Review	Package submitted; wait for the Reading Center
Changes Requested	Address Reading Center feedback and resubmit
Certified	Ready (with site Active) for production work
Not Applicable	Skip this chapter for your study



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	16 of 50

Title: NORC Portal User Guide




Figure 11. Site certification tab for studies that don't require site certification



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	17 of 50

Title: NORC Portal User Guide




Figure 12. Site certification tab for studies that require site certification

4.2.2 Recommended Certification Pathway

A typical pathway for navigating certification in the portal is as follows:

- 1. Review the Imaging Acquisition Manual (from Documents / study materials)
- 2. Enter **Site & device** information
- 3. Complete **Test upload** slots
- 4. Complete **Attestation** and submit
- 5. Wait for **Reading Center review**
- 6. Reach **Certified** (and site **Active**)

You may also receive in-app notifications when certification is approved or when changes are requested. Check the notification bell for details and reviewer comments.

4.2.3 Site and Device Information

Steps:

- 1. Go to **Certification** and open the **Site information** / site & device step. (Figure 13)



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	18 of 50

Title: NORC Portal User Guide

2. Confirm or enter the principal investigator (PI), coordinator, and related site contacts as prompted.
3. Add imaging **devices** and modalities used at your site.
4. Save and continue.

Figure 13. Site & Device details step on site certification wizard

4.2.4 Test Upload

Purpose: Upload non-production sample images into the designated certification slots to verify image quality and system functionality. These images are for certification purposes only and must not contain any subject data.

Steps:

1. Open **Certification** → **Test upload**. (Figure 14)
2. For each required slot, select and upload the sample file.
3. Wait until each required slot shows as uploaded successfully.
4. Continue to **Attestation**.



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	19 of 50

Title: NORC Portal User Guide

Figure 14. Test Upload step on site certification wizard

Notes:

- ☐ Complete all required slots before attestation.
- ☐ Use files that match the modality guidance in the Imaging Acquisition Manual.

4.2.5 Attestation and Submission of Images

Purpose: Review and confirm that all certification requirements have been completed before submitting the site certification package for reading center review.

Steps:

1. Open **Certification** → **Attestation**. (Figure 15)
2. Read and check each attestation statement.
3. Enter your name / signature fields as prompted.
4. Submit the certification package for Reading Center review.



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	20 of 50

Title: NORC Portal User Guide




Figure 15. Attestation step on site certification wizard

If status becomes **Changes Requested**, open Certification, address the feedback, and resubmit.

5. Subjects

5.1 Subjects List

Purpose: View, search, and access subjects registered at your site for the selected study.

Steps:

1. Open **Subjects**. (Figure 16)
2. Search by Subject ID or screening number as needed.
3. Filter by enrollment status if available.
4. Select a subject row to open details.



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	21 of 50

Title: NORC Portal User Guide




Figure 16. Subjects Dashboard page

5.2 Register a Subject

Purpose: Create a new subject record. Registration requires the site to be **Active** and certification complete for the specific protocol (unless certification is **Not Applicable**).

Steps:

- 1. Open **Subjects** → **Register Subject** (or use the Dashboard action). (Figure 17)
- 2. Choose how the Subject ID is assigned:
 - **Automatic** — typically **StudyCode-SiteNumber-Sequence**
 - **Protocol-assigned** — enter the ID provided by the study
- 3. Enter the screening number, screening date, and other required fields.
- 4. Set **enrollment status** and **study eye** as applicable.
- 5. Select **Register Subject**.



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	22 of 50

Title: NORC Portal User Guide




Figure 17. Register Subject screen

5.3 Subject Details and Status Updates

Purpose: View the subject visit schedule, initiate image uploads, and update enrollment when the subject’s status changes in the study.

Enrollment statuses you will see:

Status	Typical meaning
Screening	Screening period; only Screening visit imaging until Enrolled
Enrolled	Protocol follow-up visits available per schedule
Screen Failed	Terminal — no further production imaging
Withdrawn	Terminal — no further production imaging
Completed	Terminal — no further production imaging

Steps (update enrollment):



MANUAL

Document Number: MAN-018-01

Effective Date: 23-AUG-2026

Version Number: 02

Page: 23 of 50

Title: NORC Portal User Guide

1. Open the subject.
2. Change enrollment status as required by protocol. (Figure 18)
3. Confirm with your password when prompted.
4. Save.




Figure 18. Change subject enrollment status window on subject dashboard

Important: While status is **Screening**, only the **Screening** visit can be uploaded. Set status to **Enrolled** before scheduled follow-up imaging.

5.4 Subject Visits

From subject details you can typically: (Figure 19)

Action	When to use
Start Visit / Upload	Begin imaging for an available protocol visit
View Submission	Open an existing submission for that visit
Continue Upload	Finish a visit that is Partially Submitted (when the study allows partial submissions)



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	24 of 50

Title: NORC Portal User Guide

Action	When to use
Add Manual Visit	Extra visit not on the fixed schedule (as per study design)
Add Unscheduled Visit	Unscheduled imaging when the study allows it (as per study design)

Test301-101-0002

Screening SCR-0002

Enrollment status

EnrolledChange

Study eye

Right eye (OD)

Date registered

Aug 5, 2026

Open queries

0

Visit schedule

Add Manual Visit

Add Unscheduled Visit

Visit	EXPECTED	UPLOAD WINDOW	ACTUAL VISIT	DEADLINE	UPLOAD STATUS	QC	QUERIES	ACTION
Screening	Aug 5, 2026	Aug 5, 2026 – Aug 5, 2026	Aug 5, 2026	Submitted On Time	Accepted	Passed	0	View Submission
Month 4	Dec 3, 2026	Nov 30, 2026 – Dec 6, 2026	—	114 days remaining	Not Due	Not Started	0	
Month 8	Apr 2, 2027	Mar 30, 2027 – Apr 5, 2027	—	234 days remaining	Not Due	Not Started	0	
Month 12	Aug 5, 2027	Aug 2, 2027 – Aug 8, 2027	—	359 days remaining	Not Due	Not Started	0	
Unscheduled	—	—	—	Available	Available	Not Started	0	Start Visit

Figure 19. Visit Schedule for a registered subject on subject dashboard

Visit timing cues: Visits may show as **Not Due**, **Available**, **Due Soon**, or **Overdue**. Overdue uploads may require a late reason during submission.

6. Imaging submission

You can start imaging in two ways:

- ☐ **New Submission** in the sidebar (wizard)
- ☐ **Start Visit / Upload** from the subject's visit row

Both paths lead to file upload, review and attestation, and a confirmation number.

6.1 Prerequisites checklist

Before starting verify the following:

- ☐ Site status is **Active**
- ☐ Certification is **Certified** or **Not Applicable**
- ☐ Required documents are **Acknowledged** (if study policy is Required)



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	25 of 50

Title: NORC Portal User Guide

- ☐ Subject exists and enrollment allows the visit (Screening vs Enrolled)
- ☐ Files match protocol slots (modality / eye) per the Imaging Acquisition Manual

6.2 New Submission Wizard

Purpose: Upload a new imaging visit for a subject or replace existing visit imaging when allowed as per study design.

Steps:

1. Open **New Submission**.
2. Choose **submission type**:
 - Screening visit
 - Scheduled follow-up visit
 - Unscheduled visit (**per study design**)
 - Replacement images for an existing submission (**Per study design**)
3. Select the **subject**.
4. Select the **visit**.
5. Enter the **visit date**. If the visit is overdue, provide the **late reason** when prompted. (Figure 20)
6. Continue to imaging upload.

Figure 20. Visit date step on new imaging submission (start visit on subject dashboard)

Notes:

- ☐ Screening subjects must use **Screening visit** (or update enrollment first).



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	26 of 50

Title: NORC Portal User Guide

- Manual visits, when allowed, are usually started from subject dashboard, not from new submission wizard listed above.
- Links from the Dashboard attention list can open New Submission with the subject and visit already selected.

6.3 Upload Imaging Files

Purpose: Attach files to each required (and optional) imaging slot for the visit.

Steps:

- On the imaging upload page, review the list of slots (modality, eye, required or optional).
- For each required slot, choose the file and upload it. (Figure 21)
- For manual-style uploads (when applicable), enter eye, modality, and visit date as required.
- Resolve any upload errors before continuing.
- Proceed to **Review**.




Figure 21. Upload imaging step on imaging submission wizard

Supported file types (typical): .dcm, .E2E, .zip, .jpg, .jpeg, .png, .tif, .tiff — confirm any study-specific restrictions in the Imaging Acquisition Manual.



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	27 of 50

Title: NORC Portal User Guide

6.4 Partial Submissions**

If permitted by the study, a visit may be submitted when required imaging is incomplete. Enter a reason for the missing imaging, complete the attestation, and submit. (Figure 22) The visit will display as Partially Submitted. Use Continue Upload to add the remaining required images later.

***Dependent per protocol*

Figure 22. Reason for missing required imaging when submitting partial submissions on imaging upload step

After a Partial Submission:

- ☐ The visit status changes to **Partially Submitted**, and the submission is received by the Reading Center.
- ☐ To add remaining required images, select **Continue Upload** from the subject visit table or **Finish ... submission** from the Dashboard.
- ☐ Previously uploaded files are retained. Once all required imaging is submitted, the partial submission is superseded.

6.5 Review, Attest, and Submit Imaging

Steps:

1. Review validation results (Passed / Warning / Failed). Fix **Failed** items before you submit. (Figure 23) Warnings may be acceptable when the study allows partial submission, and a missing-imaging reason has been provided.
2. Read and accept the attestation statements.
3. Enter your password to confirm.
4. Submit to the Reading Center.
5. Record the **confirmation number** on the confirmation page. (Figure 24)



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	28 of 50

Title: NORC Portal User Guide




Figure 23. Upload validation and review step on imaging submission wizard



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	29 of 50

Title: NORC Portal User Guide




Figure 24. Imaging submission confirmation screen

After submission: Track the status under **Upload History** and on the subject dashboard. You should also receive an email acknowledgment that the submission was received. The Reading Center performs technical QC; you may receive queries if additional action is needed. QC outcomes (for example, accepted or changes requested) typically appear in the in-app notification bell as well as on the submission record.

6.6 Replacement Images

Use submission type **Replacement images for an existing submission**, or respond to a query with **Upload replacement image**, when the Reading Center requests re-acquisition. Follow the same upload → review → submit pattern and keep the confirmation number and query linkage for your site records.

6.7. Upload History and Submission Details

Review prior submissions, check status, and access any Reading Center correspondence.

Steps:

1. Open **Upload History**. (Figure 25)
2. Locate the submission by subject, visit, date, or confirmation number.
3. Open the submission to view status, QC information, timeline or events, and related queries. (Figure 26)



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	30 of 50

Title: NORC Portal User Guide




Figure 25. Upload history tab showing all uploads from the site




Figure 26. View submission page

Statuses you may see: Draft, Submitted, Under Technical QC, Accepted, Changes Requested, Partially Submitted (on the visit), and others listed in [Appendix A](#).

7. Automatic Slot Check

Purpose: Slot Check confirms that the file you attached **matches** the imaging slot requirements—for example the correct eye (laterality), modality, resolution, and that the file contains only one scan. It runs as soon as you



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	31 of 50

Title: NORC Portal User Guide

select a file on an imaging slot where Slot Check is **enabled**. It does not replace Automatic Grade Check, technical QC, queries, or grader review.

Slot Check is **as per study design**. When it is enabled for a slot, use the statuses and messages on the imaging upload page to correct mismatches while the patient is still available when possible.

When both Slot Check and Grade Check are enabled for the same slot, **Slot Check runs first**. Grade Check starts only after Slot Check shows **Passed**.

7.1 What you will see on imaging upload

After you select a file on a slot with Slot Check enabled:

- **Uploaded** means the file is attached to that slot.
- **Slot Check** is a separate label with its own status next to the file.

Slot Check status	Meaning	Site Action
Processing	The automatic check is still running. A spinner appears on the status. Slot Check is usually quick.	Stay on the imaging upload page briefly if you can. You do not have to wait forever—you may continue (see below).
Passed	The file matched the slot requirements (eye, modality, resolution, and single-scan rules as configured).	No action. If Grade Check is also enabled for the slot, it starts after Slot Check passes.
Failed	The file did not match one or more slot requirements. The portal shows the failure reason under the slot.	Fix the issue (see the reason table below) and use Re-upload on that slot.
Failed to process	The automatic slot check could not be completed on this file.	Verify the file, re-export if needed, then use Re-upload. You can continue to submit even if the slot check has failed to process. Additionally, you can contact support.

Re-upload appears on a slot when Slot Check is **Failed** or **Failed to process**. Selecting a new file starts Slot Check again for that slot.

There is **no separate Slot Check sidebar tab**. Review Slot Check only on the imaging upload page for that visit.



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	32 of 50

Title: NORC Portal User Guide

7.2 Failure reasons you may see

When Slot Check is **Failed**, the message under the slot is the exact reason returned by the slot check. Use this table to interpret it and correct the upload:

If you see this message	What it means	Site Action
only one scan per slot is allowed	The uploaded file contains more than one scan. Each slot expects a single scan.	Export or select a file that contains only one scan for this slot, then Re-upload .
slot modality and scan modality doesn't match	The scan's modality (for example OCT vs another type) does not match the modality required for this slot.	Confirm the slot's modality in the Imaging Acquisition Manual, acquire or select the correct modality file, then Re-upload .
slot laterality and scan laterality doesn't match	The scan's eye (OD / OS / OU) does not match the eye required for this slot.	Confirm the correct eye for the slot, acquire or select the matching eye, then Re-upload .
slot resolution and scan resolution doesn't match	The scan's resolution does not match what the slot requires (for example 97 vs 49 B-scans, or another protocol-defined count).	Re-acquire using the resolution specified for that slot in the Imaging Acquisition Manual, then Re-upload .

If Slot Check is **Failed to process**, you may see a message such as that the check did not return a result in time, or that the result could not be read. You can continue to submit or Re-upload after verifying the file; contact support if the problem continues.

7.3 Submitting while Slot Check is still running or has failed

You are not blocked from **Review & Submit**.

- If any attached slot is still **Processing** (Slot Check and/or Grade Check), you will see a confirmation that Automatic Slot Check or Grade Check is still in process. You may proceed. **Prefer** waiting briefly on the imaging page for Slot Check to finish when you can.



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	33 of 50

Title: NORC Portal User Guide

- If any attached slot is **Failed** or **Failed to process**, you will see: *One or more imaging slots failed automatic slot or grade check. Do you want to proceed with submission?* Prefer **Re-upload** after correcting the file. Proceed only if you must submit now (for example, the visit window is closing).

If both apply, one confirmation covers both messages. Choose **Stay on this page** to re-upload, or **Proceed** to continue to review, attestation, and password confirmation as usual.

Notes:

- Slot Check is enabled only for slots configured that way for your study (**as per study design**). If a slot has no Slot Check label, this feature does not apply to that slot.
- Slot Check runs per imaging slot, not on the whole visit at once.
- When Slot Check and Grade Check are both enabled, fix Slot Check failures first; Grade Check will not start until Slot Check **Passed**.
- Confirm eye, modality, and resolution against the Imaging Acquisition Manual before uploading to reduce Slot Check failures.

8. Automatic Grade Check

Purpose: Automatic Grade Check is a Reading Center quality check that evaluates imaging files for basic image quality and gradability when they are uploaded. Results are provided before technical QC to help sites identify images that may need to be reacquired or re-uploaded.

Automatic Grade Check does not replace Reading Center technical QC, queries, or grader review. When possible, review Grade Check Results while the subject is still available so that failed images can be reacquired. (Figure 27)



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	34 of 50

Title: NORC Portal User Guide




Figure 27. Imaging upload step of imaging submission wizard showing a scan that failed automatic grade check with options to re-upload the imaging

8.1 What you will see on imaging upload

After you select a file on a slot:

- ☐ **Uploaded** means the file is attached to that slot.
- ☐ **Grade Check** is a separate label with its own status next to the file.

Grade Check status	Meaning	Site Action
Processing	The automatic quality check is still in progress.	You may continue with the submission without waiting for the check to complete in case the processing takes a long time. Results can be reviewed on the imaging upload page or under Grade Check in the left menu.
Passed	The file met automatic quality requirements for grading.	No action. Continue with remaining slots or proceed to Review & Submit.



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	35 of 50

Title: NORC Portal User Guide

Grade Check status	Meaning	Site Action
Failed	The file did not meet automatic quality requirements for grading	When possible, reacquire the image and select Re-upload to replace the file.
Failed to process	The automatic check could not be completed on this file.	Re-take or re-export the file if needed, then use Re-upload. If not resolved, contact support.

Re-upload appears on a slot when Grade Check is **Failed** or **Failed to process**. Selecting a replacement file automatically begins a new Grade Check for that imaging slot.

8.2 Submitting while Grade Check is still running or has failed

- ☐ Automatic Grade Check does not prevent a site from proceeding to **Review & Submit**.
- ☐ If one or more files are still **Processing**, the portal displays a confirmation message indicating that Automatic Grade Check is still in progress. You may proceed with the submission and review the final results later under **Grade Check**.
- ☐ If one or more files show **Failed** or **Failed to Process**, the portal displays a confirmation message before allowing you to proceed. When possible, select **Stay on this page** and use **Re-upload** to replace the affected file.
- ☐ If submission cannot be delayed (for example, because the visit window is closing), select **Proceed** to continue to review, attestation, and submission.
- ☐ If both conditions apply, the portal displays a single confirmation message addressing the outstanding Grade Check results.

8.3 Grade Check tab (left menu)

Purpose: See Automatic Grade Check results for your site across subjects, visits, and imaging slots, including files you already submitted while a check was still running.

Steps:

1. Open **Grade Check** in the left sidebar. (Figure 28)
2. Use the status checkboxes to show **Passed**, **Failed**, **Processing**, and/or **Failed to process**.
3. Use **Refresh** in the header when you want the latest statuses.
4. For **Failed** or **Failed to process**, select **Re-upload** to open imaging upload for that subject, visit, and slot.
5. For other statuses, select **View upload** to open the applicable imaging submission.



MANUAL

Document Number: MAN-018-01

Effective Date: 23-AUG-2026

Version Number: 02

Page: 36 of 50

Title: NORC Portal User Guide




Figure 28. Grade Check tab on the left menu showing status of grade check for all uploads

Failed Grade Check items may also appear in **What needs attention** on the Dashboard. Selecting the item opens the same imaging upload page.

Notes:

- Grade Check runs per imaging slot, not on the whole visit at once. One slot can pass while another is still processing or has failed.
- Automatic Grade Check is immediate quality control for the site. Reading Center technical QC and graded reads still occur after the submission is received.

9. Queries

Purpose: Manage and respond to Reading Center feedback, including requests for clarification or replacement images for a submission.

9.1 Queries List

Steps:

- Open **Queries**. (Figure 29)
- Use the tabs as available: **Open**, **Awaiting Reading Center**, **Closed**, **All**.
- Click Open under Action in a query row to open query details.



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	37 of 50

Title: NORC Portal User Guide




Figure 29. Queries dashboard

When the Reading Center opens a new query, the original submitter typically receives an email, and all users at the site may see an in-app notification such as **New query received**. Open queries also appear in the Dashboard **What needs attention** list.

9.2 Respond to Queries

Steps:

- 1. Open the query and read the Reading Center request (type and comments). (Figure 30)
- 2. Choose a **response category**: Provide written response, Upload replacement image, Confirm image cannot be repeated, or Request clarification
- 3. Enter your written response.
- 4. If you are uploading replacements, attach files to the indicated slots.
- 5. Enter your name, date, and **password**, then submit.



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	38 of 50

Title: NORC Portal User Guide




Figure 30. Open Query page showing query details and options to submit response

After You Respond:

The query status changes to **Awaiting Reading Center** while your response is reviewed. Once resolved, the query moves to **Closed**. If additional information is needed, the Reading Center may reopen the query.

Notes:

- ☐ Select the response category that best matches the Reading Center’s request.
- ☐ Retain confirmation numbers and query IDs in accordance with your site SOPs.

10. Data Dumps

Purpose: Transfer an unstructured .zip package to the Reading Center outside the per-visit imaging wizard. Data dumps are **per study design** — the **Data Dumps** menu appears only when the study has data dumps enabled. If you do not see the menu, this feature is not available for your study or the site may not be listed as Active.

10.1 Create and Upload a Data Dump

Steps:



MANUAL

Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	39 of 50

Title: NORC Portal User Guide

1. Open **Data Dumps**. (Figure 31)
2. Select **New** / create a new data dump. (Figure 32)
3. Choose a **.zip** file only.
4. Write a short description of data dump.
5. Confirm with your password when prompted and start the upload.
6. Remain on or return to the dump **detail** page to monitor progress. Large files may continue uploading in the background while that page is open. (Figure 33)
7. When finished, status should show **Submitted**. Keep the confirmation number if shown.




Figure 31. Data dumps dashboard




Figure 32. New data dump screen



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	40 of 50

Title: NORC Portal User Guide




Figure 33. Data dump status page

Statuses: Uploading, Paused, Failed, Submitted. If upload fails or pauses, retry from the detail page or contact support with the dump ID or confirmation number.

11. Support and Communications

11.1 Study Contacts

Purpose: **Study Contacts** provides the Reading Center and technical contacts available for your study, including names, roles, email addresses, phone numbers, and hours when provided. Select an email address to open your default email application.

If no contacts are configured for the study, you will see a message such as **No contacts configured for this study**. The Dashboard may also offer a quick view of the primary study contact and imaging manual reference.

11.2 Help (FAQs and Contact support)

Purpose: **Help** provides FAQs for common portal tasks and a **Contact Support** form for sending questions directly to the support team. Replies are sent to the email address associated with your portal account. If you cannot access the portal at all, use the offline contacts in [Appendix C](#).



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	41 of 50

Title: NORC Portal User Guide




Figure 34. Help tab on left menu

For protocol imaging technique questions, always prefer the study-specific *Imaging Acquisition Manual and Documents* over the short FAQ text.

Steps for contacting support:

1. Open **Help** and go to the **Contact support** section. (Figure 34)
2. Read the guidance: you are sending a message to the support team, and replies will go to your account email.
3. Enter your **Message** in the text area (describe your question or issue).
4. Select **Send message**. While the message is sending, the button may show **Sending....**
5. On success, you should see a confirmation such as **Help message sent**, and the form clears.

11.3. Notifications

Purpose: Alerts for queries, certification updates, QC status changes, new documents, and other site-relevant events.

Steps:

1. Select the **bell** icon in the header. (Figure 35)
2. Open a notification to go to the related page (query, certification, submission, and so on).



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	42 of 50

Title: NORC Portal User Guide

3. Optionally select **Mark all as read**.




Figure 35. Notification menu (bell icon)

Examples of in-app notifications include:

- ☐ New query received
- ☐ Query closed
- ☐ Upload accepted / submission under review / changes requested on a submission
- ☐ Site certification approved / certification changes requested
- ☐ New study document published (when acknowledgment is required)

Check notifications daily during active enrollment and imaging windows.

11.4 Email notifications



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	43 of 50

Title: NORC Portal User Guide

Purpose: The portal provides in-app and email notifications for important account and study events. Select the **bell icon** to view in-app notifications and open the related item. Monitor the email address associated with your portal account for verification codes, account updates, submission acknowledgments, and query notifications.

For a related workflow summary, see [Imaging submission](#) and [Queries](#).



Title: NORC Portal User Guide

MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	44 of 50

12. Revision History

Version Number	Date	Summary of Change(s)
01	13-AUG-2026	New Document
02	23-AUG-2026	Addition of Section 7 ‘Slot Check’



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	45 of 50

Title: NORC Portal User Guide

13. Signature Page

Author:

NAME	TITLE	SIGNATURE	DATE
Sharat Chandra Vupparaboina	Data Engineer	Digitally signed by Name : Sharat Chandra Vupparaboina Email : sharat@netramind.ai Reason : I reviewed and agree to this document Time : Aug 22 2026 20:59 EDT Signed with: Zoho Sign	Aug 22 2026

QA:

NAME	TITLE	SIGNATURE	DATE
Justin Shaka	Fractional COO	Digitally signed by Name : Justin Shaka Email : justin@netramind.ai Reason : I reviewed and agree to this document Time : Aug 22 2026 20:52 EDT Signed with: Zoho Sign	Aug 22 2026

Approver:

NAME	TITLE	SIGNATURE	DATE
Sandeep Chandra Bollepalli	Director of AI	Digitally signed by Name : Sandeep Chandra Bollepalli Email : sandeep@netramind.ai Reason : I approve this document Time : Aug 22 2026 20:55 EDT Signed with: Zoho Sign	Aug 22 2026



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	46 of 50

Title: NORC Portal User Guide

Appendix A — Key Status reference

A.1 Site Certification

Status	Meaning/Site Action Needed
Pending Review	Certification has been submitted and is awaiting Reading Center review.
Changes Requested	Reading Center feedback requires site action; review the feedback, make the requested updates, and then resubmit.
Certified	Site certification has been approved by the Reading Center.
Not Applicable	Site certification is not required for the study.

A.2 Imaging Submissions

Status	Meaning/Site Action Needed
Draft	The submission has been started but has not yet been submitted; complete and submit the imaging when ready.
Partially Submitted	The submission was sent with required imaging outstanding; use “Continue Upload” to add the remaining images.
Submitted	The submission has been received by the Reading Center.
Under Technical QC	The Reading Center is reviewing the submitted imaging.
Accepted	The submission has completed technical QC and has been accepted.
Changes Requested	Additional site action is required. Review the associated feedback or query.



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	47 of 50

Title: NORC Portal User Guide

Superseded	The submission is no longer the current version, such as when a partial submission has been completed with additional imaging.
------------	--

A.3 Queries

Status	Meaning/Site Action Needed
Open	A response or other action is required from the site.
Awaiting Reading Center	The site has responded, and the query is awaiting Reading Center review.
Closed	The query has been resolved. No further site action is required unless it is reopened.

Appendix B — Glossary

Term	Definition
NORC / Reading Center	Netramind Ophthalmic Reading Center — reviews certification, QC, and queries
NORC Site Portal	Clinical site-facing application described in this guide
MFA / verification code	One-time code emailed during login, signup, or password reset
Site certification	Site readiness package (devices, test images, attestations)
Imaging Acquisition Manual	Study imaging procedure manual (versioned)
Document acknowledgment	Site confirmation that required study documents were reviewed
Subject ID	Unique subject identifier in the study
Screening number	Site screening identifier, when used by the study
Study eye	Right / Left / Both / Not yet assigned; files may use OD / OS / OU
Imaging slot / modality	Required or optional image type for a visit (for example OCT, FAF)



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	48 of 50

Title: NORC Portal User Guide

Term	Definition
Visit window / deadline	Protocol timing for when imaging should be uploaded
Manual visit	Site-created visit outside the standard named schedule (as per study design)
Unscheduled visit	Off-schedule visit when the study allows it
Partial submission	Protocol visit package submitted with some required slots still empty, when the study allows it
Continue Upload	Action to add remaining imaging to a Partially Submitted visit
Attestation	Signed confirmations before certification or production submit
Technical QC	Reading Center quality check of submitted imaging
Query	Reading Center request for clarification or replacement images
Replacement images	Re-upload for a prior submission or query
Data dump	Unstructured .zip transfer when enabled for the study
Study Contacts	Directory of published Reading Center and support contacts for the study

Appendix C — Support and troubleshooting

C.1 Support contacts

When you can access the portal, refer to:

- **Study Contacts** for published names, roles, email (mail to), phone, and hours
- **Help → Contact support** to send a message from the portal (replies go to your account email)

If you cannot access the portal:

Channel	Contact
Email	norc@netramind.ai



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	49 of 50

Title: NORC Portal User Guide

Include your **site number**, **study code**, user email, and any **confirmation number** or **query ID** when you contact support.

C.2 Common issues

Issue	What to try
Cannot sign in after signup	Wait for account approval; confirm email; reset password only after the account is active. Watch for approval or rejection email.
MFA code not received	Check spam; wait and request a new code; confirm inbox access, contact norc@netramind.ai
Signed out unexpectedly	Check whether the login page shows an inactivity message; sign in again. Use Stay signed in if the idle warning appears.
Cannot register subjects	Confirm site is Active and certification is Certified or Not Applicable
Cannot upload imaging	Confirm document acknowledgment (if required), site Active, enrollment allows the visit, and no blocking open query on that visit
Cannot submit with missing slots	Either complete all required slots, or—if the study allows partial submission—provide a Reason for missing required imaging and submit as partial
Finish a partial visit	Use Continue Upload on the subject visit row, or the Dashboard attention item Finish ... submission
Screening follow-up blocked	Set enrollment to Enrolled before non-screening visits
Overdue visit	Provide late reason in the submission wizard
Query submit fails	Complete required text/files, password, and category
Did not receive expected email	Check spam and the inbox used for the portal account; also check the notification bell for in-app alerts
Contact support send fails	Enter at least about 10 characters; try again; if the portal is unavailable, use Study Contacts or offline support



MANUAL	
Document Number:	MAN-018-01
Effective Date:	23-AUG-2026
Version Number:	02
Page:	50 of 50

Title: NORC Portal User Guide

Issue	What to try
Data Dumps missing	Feature not enabled for this study (as per study design)
Data dump failed	Retry from dump detail; confirm .zip only; contact support with dump ID
Certification Changes Requested	Open Certification, address feedback, resubmit

C.3 Related study materials to reference

- ☐ Imaging Acquisition Manual (version listed under study documents)
- ☐ Protocol visit schedule and imaging requirements for your study